



Clarity™ ADVANCED Ceramic Brackets Parts List

CLARITY™ | ADVANCED
advanced ceramic brackets

Maxillary			IN/OUT		0.018 in. Slot		0.022 in. Slot	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Central	+17	+4	0.026	0.66	006-201	006-202	006-301	006-302
Lateral	+10	+8	0.035	0.89	006-205	006-206	006-305	006-306
Cuspid HK	0	+8	0.021	0.53	006-209	006-210	006-309	006-310
Universal Bicuspid	-7	0	0.029	0.74	006-211		006-311	
Mandibular			IN/OUT		0.018 in. Slot		0.022 in. Slot	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Anterior	-6	0	0.045	1.14	006-250		006-350	
Cuspid HK	0	+3	0.020	0.51	006-253	006-254	006-353	006-354
1st Bicuspid	-12	+2	0.030	0.76	006-255	006-256	006-355	006-356
2nd Bicuspid	-17	+2	0.034	0.86	006-259	006-260	006-359	006-360
Adding Prefix "3" is APC™ II Adhesive System								
Adding Prefix "5" is APC™ PLUS Adhesive System								

Single Patient Kits	0.018 in. Slot	0.022 in. Slot
U/L 5x5 Clarity™ ADVANCED Brackets Cuspid HK	006-100	006-110
U/L 3x3 Clarity ADVANCED Brackets Cuspid HK	006-103	006-113
U 5x5 Clarity ADVANCED Brackets Cuspid HK	006-105	006-115
U 3x3 Clarity ADVANCED Brackets Cuspid HK	006-108	006-118
Adding Prefix "3" is APC II Adhesive System		
Adding Prefix "5" is APC PLUS Adhesive System		

APC™ II
Available in APC™ II Adhesive
Coated Appliance System

APC™ PLUS
Available in APC™ PLUS Adhesive
Coated Appliance System

MBT™
Versatile+ Appliance System

3M Unitek

Clarity™ ADVANCED Ceramic Brackets Practice Marketing Tools

CLARITY™ | ADVANCED
advanced ceramic brackets

Description	Part Number
Bracket Typodont	600-236
20x Demonstration Model	600-237
Patient Brochure - Large	021-118
Patient Brochure - Small	021-119
Lenticular - Comparison With Metal Braces	012-255
Desktop Display (11x8.5 in.) - Group	014-536
Desktop Display (8.5x11 in.) - Blonde Woman	014-537
Desktop Display (8.5x11 in.) - Brunette Woman	014-538
Desktop Display (8.5x11 in.) - Dark Hair Woman	014-539
Desktop Display (8.5x11 in.) - Male	014-540
Office Poster (24x18 in.) - Group	014-541
Office Poster (18x24 in.) - Blonde Woman	014-542
Office Poster (18x24 in.) - Brunette Woman	014-543
Office Poster (18x24 in.) - Dark Hair Woman	014-544
Office Poster (18x24 in.) - Male	014-546

For product and model images to use in your practice marketing, contact 3M Unitek at 800-423-4588 and ask for extension 4303.

Desktop Displays/Office Posters



3M

3M Unitek

Orthodontic Products

2724 South Peck Road
Monrovia, CA 91016 USA
www.3MUnitek.com

In U.S. and Puerto Rico: 1-800-423-4588 • 626-574-4000

In Canada: 1-800-443-1661

Technical Helpline: 1-800-265-1943 • 626-574-4577

CE Hotline: 1-800-852-1990 x4649 • 626-574-4649

Outside these areas, contact your local representative.

APC, Clarity, and MBT are trademarks of 3M.

Please recycle. Printed in USA.

© 2012 3M. All rights reserved.

012-261-1 1204