

# Buccal Tubes System



Buccal tubes with comfortable and strong construction design. These tubes are manufactured with the most resistant steel available in the market, providing reliable and predictable orthodontic solutions for safe orthodontics.

## Characteristics and Benefits

80-micron base mesh:  
Better bond strength

17/4 steel:  
Better strength

Integral hooks (one piece construction):  
No deformation and breakage

Vertical demarcation:  
Assists in accurate bracket positioning

Rounded corners:  
Enhance comfort for patients





Note: Prewelded tubes are available.  
For further information, please contact our customer service department.



Buccal Tubes welding  
Buccal Tubes bonding






Buccal Tubes - welding

U = UPPER

Tooth	Prescription	Torque	Angulation	Offset		Occlusal Gingival	Convertible	Length (mm)	Width Distal	Hdgr Tube	Slot .018		Slot .022	
				Distal							Right	Left	Right	Left
First Molar Single 	Roth	-10°	0°	10°		-	Yes	3.5	Yes	-	336-877	336-878	336-977	336-978
	MBT	-14°	0°	10°		-	Yes	3.5	Yes	-	336-873	336-874	336-973	336-974
	Edgewise - U/L - D/E	0°	0°	0°		-	Yes	4.0	Yes	-	300-201	300-201	300-301	300-301
First Molar Convertible Double 	Roth	-10°	0°	10°	0		Yes	3.5	Yes	.045	326-825	326-826	326-925	326-926
	Roth	-10°	0°	10°	G		Yes	3.5	Yes	.045	326-829	326-830	326-929	326-930
	MBT	-14°	0°	10°	0		Yes	3.5	Yes	.045	326-871	326-872	326-971	326-972
	MBT	-14°	0°	10°	G		Yes	3.5	Yes	.045	326-873	326-874	326-973	326-974
	Rickets - U/L - D/E	0°	0°	0°	O/G		No	3.5	Yes	-	342-011	342-011	-	-
	Edgewise	0°	0°	0°	0		Yes	3.5	Yes	.045	326-805	326-806	326-905	326-906
First Molar Convertible Triple 	Roth	-10°	0°	10°	0		Yes	3.5	Yes	.045	340-821	340-822	340-921	340-922
	MBT	-14°	0°	10°	0		Yes	3.5	Yes	.045	340-873	340-874	340-973	340-974
	Rickets	0°	0°	0°	G		Yes	3.5	Yes	.045	340-011	340-012	-	-
Second Molar Single 	Roth	-10°	0°	10°	-		No	3.5	Yes	-	324-211	324-212	324-311	324-312
	MBT	-14°	0°	10°	-		No	3.5	Yes	-	324-273	324-274	324-373	324-374
	Edgewise	0°	0°	0°	-		No	3.5	Yes	-	324-001	324-002	324-101	324-102





# Buccal Tubes System

L = LOWER

Tooth	Prescription	Torque	Angulation	Offset		Occlusal	Convertible	Length (mm)	Width Distal	Hdgr Tube	Slot .018		Slot .022	
				Distal	Gingival						Right	Left	Right	Left
First Molar Single 	Roth	-25°	0°	5°	-	Yes	3.5	Yes	-		336-813	336-814	336-913	336-914
	MBT	-20°	0°	0°	-	Yes	3.5	Yes	-		336-815	336-816	336-915	336-916
First Molar Convertible Double   	Roth - Retangular	-25°	0°	5°	G	Yes	3.5	Yes	-		342-813	342-814	342-913	342-914
	Roth - Pla	-25°	0°	5°	G	Yes	3.5	Yes	.045		339-813	339-814	339-913	339-914
	MBT - Retangular	-20°	0°	0°	G	Yes	3.5	Yes	-		342-815	342-816	342-915	342-916
	MBT - Pla	-20°	0°	0°	G	Yes	3.5	Yes	.045		339-815	339-816	339-915	339-916
	Rickets - U/L - ED	0°	0°	0°	O/G	Não	3.5	Yes	-		342-011	342-011	-	-
Second Molar Single 	Roth	-25°	0°	5°	-	No	3.5	Yes	-		324-231	324-232	324-331	324-332
	MBT	-10°	0°	0°	-	No	3.5	Yes	-		324-293	324-294	324-393	324-394
	Edgewise	0°	0°	0°	-	No	3.5	Yes	-		324-001	324-002	324-101	324-102






Buccal Tubes bonding

U = UPPER

Tooth	Prescription	Torque	Angulation	Offset		Occlusal	Convertible	Length (mm)	Width	Hdgr Tube	Slot .018		Slot .022	
				Distal	Gingival						Right	Left	Right	Left
First Molar Single 	Roth	-10°	0°	10°	-	-	Yes	3.5	Yes	-	351-821	351-822	351-921	351-922
	MBT	-14°	0°	10°	-	-	Yes	3.5	Yes	-	351-873	351-874	351-973	351-974
	Edgewise - U/L - D/E	0°	0°	0°	-	-	Yes	4.0	Yes	-	301-201	301-201	301-301	301-301
First Molar Convertible Double 	Roth	-10°	0°	10°	0	0	Yes	3.5	Yes	.045	351-832	351-833	351-932	351-933
	Roth	-10°	0°	10°	G	0	Yes	3.5	Yes	.045	351-823	351-824	351-930	351-931
	MBT	-14°	0°	10°	0	0	Yes	3.5	Yes	.045	351-834	351-835	351-934	351-935
	MBT	-14°	0°	10°	G	0	Yes	3.5	Yes	.045	351-275	351-276	351-375	351-376
	Rickets - U/L - D/E	0°	0°	0°	0/G	0	No	3.5	Yes	-	351-015	351-015	-	-
	Edgewise	0°	0°	0°	0	0	Yes	3.5	Yes	.045	351-805	351-806	351-905	351-906
First Molar Convertible Triple 	Roth	-10°	0°	10°	0	0	Yes	3.5	Yes	.045	351-826	351-827	351-926	351-927
	MBT	-14°	0°	10°	0	0	Yes	3.5	Yes	.045	351-877	351-878	351-920	351-919
	Rickets	0°	0°	0°	G	0	Yes	3.5	Yes	.045	351-011	351-012	-	-
Second Molar Single 	Roth	-10°	0°	10°	-	-	No	3.5	Yes	-	351-633	351-634	351-733	351-734
	MBT	-14°	0°	10°	-	-	No	3.5	Yes	-	351-273	351-274	351-373	351-374
	Edgewise	0°	0°	0°	-	-	No	3.5	Yes	-	351-001	351-002	351-101	351-102

# Buccal Tubes System

L = LOWER

Tooth	Prescription	Torque	Angulation	Offset		Occlusal	Convertible	Length (mm)	Width Distal	Hdgr Tube	Slot .018		Slot .022	
				Distal	Gingival						Right	Left	Right	Left
First Molar Single 	Roth	-25°	0°	5°	-	Yes	3.5	Yes	-		351-813	351-814	351-913	351-914
	MBT	-20°	0°	0°	-	Yes	3.5	Yes	-		351-815	351-816	351-915	351-916
First Molar Convertible Double   	Roth - Retangular	-25°	0°	5°	G	Yes	3.5	Yes	-		351-828	351-829	351-928	351-929
	Roth - Pla	-25°	0°	5°	G	Yes	3.5	Yes	.045		351-817	351-818	351-917	351-918
	MBT - Retangular	-20°	0°	0°	G	Yes	3.5	Yes	-		351-819	351-820	351-923	351-924
	MBT - Pla	-20°	0°	0°	G	Yes	3.5	Yes	.045		351-291	351-292	351-391	351-392
	Rickets - U/L - ED	0°	0°	0°	O/G	No	3.5	Yes	-		351-015	351-015	-	-
Second Molar Single 	Roth	-25°	0°	5°	-	No	3.5	Yes	-		351-653	351-654	351-753	351-754
	MBT	-10°	0°	0°	-	No	3.5	Yes	-		351-293	351-294	351-393	351-394
	Edgewise	0°	0°	0°	-	No	3.5	Yes	-		351-001	351-002	351-101	351-102