



American Crafted

dentronix 
Instrumentation. Sterilization. Innovation.



“

I regarded our progress merely as an invitation to do more - as an indication that we had reached a place where we might begin to perform a real service.

















Henry Ford

”

American Crafted



Dentronix is committed to homegrown products. We have been building orthodontic pliers in the USA since the late 1960's and were one of the pioneers of inserted-tip design. As inventors of many useful instrument designs, Dentronix became a leader in the infection control revolution of the 1980's, designing and developing the most popular sterilizer for orthodontists. Today, we continue our tradition of made-in-the-USA products that are manufactured to meet your needs and exceed your expectations.

instrument selection guide	2	
distal end cutters	3	
cutters	4	
wire bending	5	
utility pliers	10	
band & bracket removal	12	
band & bracket placement	13	
laboratory	14	
scalers, pushers, pluggers, explorers	16	
hemostats	17	
mirrors & scissors	18	
specialty instruments & accessories	19	
records	20	
sterilization simplified	21	
sterilizer accessories	26	
disinfectant cleaner	29	

Classic

The original-design **dentronix** pliers that made us a leader in the orthodontic instrument business. Inserted with vanadium cobalt tool steel, these pliers offer the optimum in wear, control, and durability. Three year warranty on manufacturing defects. Lifetime warranty on bond failure. Dry Heat sterilize only.

Royal

Our autoclavable line of pliers feature super-hard non-ferrous (mostly carbide) inserts that resist wear and corrosion even when steam sterilized. Lifetime warranty on bond failure and corrosion. Five year warranty against manufacturing defects. One year warranty on reconditioning, sharpening, and reinserting.

NI non-inserted

Our non-inserted line offers a steam sterilizable option made from solid stainless steel. These pliers are hardened to optimum levels and offer good wear characteristics, however not to the level of our inserted Classic and Royal lines. Non-inserted pliers have a 90-day warranty against manufacturing defects.

Comparison of Tip Materials

Properties	TOOL	Carbide	Stainless Steel
	STEEL		
Hardness	✱	✱	
Ductility	✱		✱
Edge Holding	✱	✱	✱
Shock Resistance	✱		✱
Grip	✱		✱
Red Hardness	✱	✱	
Thermal Expansion	✱		
Grindability	✱		✱
Corrosion Resistance		✱	✱


Box Joint VS. Lap Joint

All dentronix pliers are manufactured from premium USA-made box-jointed plier blanks. Compared to the scissors-style lap joint, this encapsulated design offers tremendous advantages in increased bearing surfaces and contact areas for maintaining tip alignment.



High Gloss VS. Satin Finish

Standard finish on dentronix pliers is high gloss. This highly polished surface provides the best resistance to corrosion, staining, and discoloration that is available to the orthodontic professional. We believe that for a premium plier, it is the only sensible finish. For those who prefer a low gloss satin finish, it can be provided on request at no additional charge.



DISTAL END CUTTERS

Original Shear-Hold® Mini

Classic D200
Royal DE200
NI
.022" x .025" MAX



Classic D200S
Royal DE200S
NI
.022" x .025" MAX



Original Shear-Hold®

*dentronix invented
and patented the
Shear Hold®
safety hold distal
end cutter in 1971.*

Safety Hold Flush Cut V Design

Classic D200V
Royal DE200V
NI
.018" x .022" MAX
.016" MIN



Safety Hold Flush Cut V Design Mini

Classic D200SV
Royal DE200SV
NI
.018" x .022" MAX
.016" MIN



Safety Hold Flush Cut Spring Design Mini NiTi ONLY

Classic
Royal DE200SNT
NI
.010" to .018" NiTi



Safety Hold Flush Cut Spring Design

Classic D200F
Royal
NI
.018" x .022" MAX



Non-Hold Flush Cut

Classic D199
Royal DE199
NI
.018" x .022" MAX



to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

CUTTERS



Classic D201
Royal DE201
NI
.016" Soft MAX

Ligature Cutter
Straight



Classic D201S
Royal DE201S
NI
.016" Soft MAX

Ligature Cutter
Straight Mini

Ligature Cutter 15°



Classic D201A
Royal DE201A
NI
.016" Soft MAX

Ligature Cutter 15°
Mini



Classic D201SA
Royal DE201SA
NI
.016" Soft MAX

Ligature Cutter
Tapered Mini



Classic D201ST
Royal DE201ST
NI
.016" Soft MAX

Ligature Cutter 45°



Classic D212
Royal DE212
NI
.014" Soft MAX

Hard Wire Cutter
Straight



Classic D201H
Royal DE201H
NI
.022" x .028" MAX



Serrated Bird Beak

Bird Beak



Bird Beak
with Cutter



Bird Beak
Large Tips



Light Wire
Magness



Jarabak



to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

WIRE BENDING



Light Wire with Cutter
No Grooves

Light Wire
No Grooves



Light Wire with Cutter
Two Grooves



Light Wire
Two Grooves



Light Wire with Cutter
Three Grooves



Light Wire
Three Grooves



Light Wire & Arch Contouring Plier No Grooves



Classic D209TR
Royal
NI
.022" MAX



Classic D104P
Royal DE104P
NI
.020" MAX
Tips Measure: ● — ■
.025" .025"

Roth Light Wire Design by Dr. Ron Roth

Nance Loop Plier



Classic D203
Royal DE203
NI
.022" MAX

Samson Plier



Classic D208S
Royal
NI
.028" MAX

Tweed Arch



Classic D210
Royal DE210
NI
.022" x .028" MAX

Tweed Arch Short Jaw



Classic D210T
Royal DE210T
NI
.022" x .028" MAX

Tweed Omega Loop



Classic D236
Royal DE236
NI
.022" x .028" MAX
● — ■
.045" .100"

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

WIRE BENDING

Optical Plier

Classic D208
Royal DE208
NI

.028" MAX

→ ● ← ↺
.035" .100"



Classic D242
Royal DE242
NI

.022" MAX



V Bend Stop Plier



Keyless Torquing Set

Female

Male

V Bend Stop Plier

NiTi

NiTi Wire Bender

Lingual Arch Plier

Classic D261A
Royal DE261A
NI
.022 x .025" MAX



Classic D261B
Royal DE261B
NI
.022" x .025" MAX



Classic D242N
Royal DE242N
NI
.022" MAX



Classic D222NT
Royal NI
.016" MAX



Classic D243
Royal DE243
NI
.036" MAX



WIRE BENDING

Turbyfill Loop
Impingement Plier



Classic D131
Royal
NI
.022" MAX

Detailing Plier 1mm



Classic D270-1
Royal DE270-1
NI
.022 x .028" " MAX

Detailing Plier 3/4mm



Classic D270-.7
Royal DE270-.7
NI
.022 x .028" " MAX

Detailing Plier 1/2mm



Classic D270-.5
Royal DE270-.5
NI
.022" x .028" MAX

Arch Contouring Plier
Grooved



Classic DE211
Royal
NI
.025" MAX

Arch Contouring Plier
Non-Grooved



Classic DE211P
Royal
NI
.025" MAX

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

UTILITY PLIERS



Weingart

Classic D213I
Royal DE213I
NI D213
.022" x .028" MAX

Small Weingart



Classic
Royal
NI D213S
.022 x .028" " MAX

How 3mm tips



Classic D110I
Royal
NI .022" x .028" MAX

How 4mm tips



Classic
Royal
NI D110
.022" x .028" MAX

How Angled
6mm tips



Classic
Royal
NI D3195
.022" x .028" MAX

Inserting Plier
Smooth Jaws



Classic D241
Royal DE241
NI

UTILITY PLIERS

Quad Helix Plier



Classic D344Q
Royal
NI
.038" MAX

Adjusts like a 3-prong plier without removing wire.

Lingual Arch Placer



Classic D260
Royal
NI
.036" MAX

Cinch Back Plier 90°



Classic D502
Royal DE502
NI
.022" x .028" MAX

Cinch back Plier 45°



Classic D503
Royal DE503
NI
.022" x .028" MAX

Crimpable Hook Plier
Straight



Classic D172
Royal
NI

Crimpable Hook Plier
Angled



Classic D172A
Royal
NI

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

BAND & BRACKET REMOVAL



Classic D347SI
Royal DE347SI
NI

Posterior Band Remover
Short



Classic D231
Royal DE231
NI



Anterior Debonder
& Scaler



Diamond Sharpening File for Bracket Removers
DDS-1

Posterior Band Remover
Long



Classic D347I
Royal D347
NI

Angled Debonder
Narrow



Classic D230N
Royal DE230N
NI



Angled Debonder
Wide



Classic D230W
Royal
NI

Anterior Debanding Plier



Classic D160
Royal DE160
NI

Band Slitter



Classic D215
Royal
NI

Bracket Placing Tweezers etc.

D219

D220

D226

D221

D227

D228

D229

D527

Classic
Royal
NI

D3125

High Heat Bite Sticks

DBS4
Up to
450°F

DBS2
Up to
270°F

Boone Positioning Gauge

DPG-1

Band Pusher

D300

Dickenson Band Contouring Plier

BAND & BRACKET PLACEMENT

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

Bird Beak Standard



Classic D139
Royal DE139
NI

.030" MAX

Tips Measure: ● ■
.035" .035"



Classic D134
Royal DE134
NI

.030" MAX

Tips Measure: ● ■
.028" .035"

Serrated Bird Beak

Three Prong Medium Tips



Classic D344
Royal
NI .036" MAX

Three Prong Fine Tips



Classic D344T
Royal
NI .030" MAX

Large Bird Beak Medium Tips



Classic D134M
Royal DE134M
NI .036" MAX
Tips Measure: ● ■
.032" .040"

Large Bird Beak Large Tips



Classic D134L
Royal DE134L
NI .056" MAX
Tips Measure: ● ■
.045" .055"

Young Plier



Classic D3144
Royal
NI .045" MAX



Facebow Plier Set

D280SET

Wire Cutter

D281

.070" MAX

Wire Bender/ Former

D282

Compound Action Cutter

D177

.051" MAX



RES

Ratcheting Expansion Screw

Innovative design prevents back-up.
15 mm width/11 mm expansion. Use
with special safety key



Bernard Plier 5-1/2"

D40A



Bernard Plier 6-1/2"

D40B



Heavy Three Prong

Classic
Royal
NI

D344F
.072" MAX



Standard Adams Plier

Classic D135
Royal DE135
NI

.030" MAX
Tips Measure: .035" .035"



Large Adams Plier

Classic D135L
Royal DE135L
NI

.051" MAX
Tips Measure: .028" .040"



LABORATORY

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

15

SCALERS PUSHERS PLUGGERS EXPLORERS



D15	Pusher/Scaler	D306	Pusher/Director	D305	Ligature Director
D15A	Plugger/Scaler	D307	#3 Scaler/Director	D307B	Scaler
D308	B Scaler/Director	D16	Pusher/Pusher	D6	Explorer
D3-3S	Scaler/Whiteside	D1612	B Scaler/45° Pusher	D23	Explorer/Sickle


Roth Band Driver



DP1001-01TR


Replacement Band Driver Tip
DP1002-01TR

Mathieu Style



Non-inserted
Tapered Tip

D216



Hole-in-Tip design

D216O



Extra Tapered Tip

D216T



Carbide Inserts
Tapered Tip

D216I



Hook Tip

D223

Mosquito Hemostats



D32H
4-1/2"
Straight Tip



D33H
4-1/2"
Curved Tip



D32HO
4-1/2"
Hole-in-Tip



D224
4-1/2"
Hook Tip

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

MIRRORS SCISSORS




Double-sided
Mirror Only
D5D



#4 Rhodium
Single-sided
D4



#5 Rhodium
Single-sided
D5



Handle Only
D45H

Double-sided
D5DC



Crown & Collar
Curved Serrated
D428



Festooning
D677

Crown & Collar
Straight
D825SS



Crown & Collar
Curved
D827SS






Cotton Plier
D317



Cam-Lok
Elastic Positioner
D225



Twirl-On
Elastic Positioner
D31



Hi-Heat Cheek Retractors
DCH1 STANDARD
DCH2 PEDO
Dry Heat Sterilizable




Separating Plier
D151

SPECIALTY INSTRUMENTS & ACCESSORIES

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com



Cephalometric Protractor
DCP-1



Consultation Typodont
DTYP

Typodont Training Aid

Frame	DT1
Metal Teeth	DT2
Normal Occlusion	DT3
Class I - Crowded	DT31
Class I - Division 1	DT321
Class II - Division 2	DT322
Class III	DT33



Disposable Impression Trays
Packs of 12

Upper/Large	DC6901
Lower/Large	DC6902
Upper/Medium	DC6903
Lower/Medium	DC6904
Upper/Small	DC6905
Lower/Small	DC6906



Upright Rack 18"
DPR-18



Designer Rack
DPR15



Upright Rack 5"
DPR-5

4" Diameter Round Rack
with center insert
DPR-1



4" Diameter
Round Rack
DPR-1A



Mixing Bowls

Small Bowl	3.25" diameter	D563
Medium Bowl	4.5" diameter	D564
Large Bowl	4.75" diameter	D565
X-Large Bowl	5" diameter	D566

Mixing Spatula
DSP1



to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com



Dentronix helps make your instrument handling process simple and efficient with products designed to work together to safeguard your investment, your staff, and your patients. All products are specified to be consistent with state and federal regulations and guidelines to afford you the highest measure of compliance.



DDUS60 Ultrasonic Cleaning System

An essential element of the dentronix dry heat sterilization process, the DDUS60 series of ultrasonic cleaners is designed to simultaneously clean four fully-loaded DDS7000 plier racks and provide a convenient hands-free means of draining excess liquid prior to loading in the sterilizer.

DDUS60 Countertop Unit Outside Dimensions

19"W x 12-1/4"D x 14-3/8" H
(48.3 cm x 31.1cm x 37.2cm)
Height with lid open - 23-1/2"

Tank Dimensions

14" W x 9"D x 8-3/4"H
(35.5cm x 22.8cm x 22.2cm)

Tank Capacity

4 Gallons

Operating Frequency

40 kHz

Shipping Weight

45 lbs

Power Output

400 Watts

All DDUS60 units include:

Hinged Cover
Timer
Front-mounted drain
Choice of Plier Rack or Deep
Basket
Two year Warranty



**Move instruments to cleaning area
decontaminate treatment area and replace barriers.**

"If manual cleaning is not performed immediately, placing instruments in a puncture-resistant container and soaking them with detergent, a disinfectant/detergent, or an enzymatic cleaner will prevent drying of patient material and make cleaning easier and less time-consuming." CDC Guidelines

DDUS60115	115V	unit with Plier Rack Stand
DDUS60115B	115V	unit with Deep Basket
DDUS60230	230V	unit with Plier Rack Stand
DDUS60230B	230V	unit with Deep Basket

2

Clean Instruments

by hand or with ultrasonic cleaner.

"Reusable instruments, supplies, and equipment should be received, sorted, cleaned, and decontaminated in one section of the processing area. Cleaning should precede all disinfection and sterilization processes; it should involve removal of debris as well as organic and inorganic contamination. Removal of debris and contamination is achieved either by scrubbing with a surfactant, detergent, and water, or by an automated process (e.g., ultrasonic cleaner or washer-disinfector) using chemical agents." CDC Guidelines



DDUS60R Recessed Ultrasonic Cleaning System



Recessed Unit DDUS60R Above Counter
19" W x 12-1/4" D x 2-1/4" H
(48.3 cm x 31.1 cm x 5.7 cm)
Height with lid open - 11.5"

Below Counter
14-5/8" x 10-5/8" x 18-3/8"
(37.1 cm x 26.9 cm x 49.3 cm)

Generator Dimensions
17" x 11-1/2" x 4-3/4"
(43.2 cm x 29.2 cm x 12.1 cm)

Tank Dimensions
14" W x 9" D x 8-3/4" H
(35.5 cm x 22.8 cm x 22.2 cm)

Shipping Weight
70 lbs

Power Output
400 Watts

DDUS60R115 115V unit with Plier Rack Stand
DDUS60R115B 115V unit with Deep Basket

DDUS60R230 230V unit with Plier Rack Stand
DDUS60R230B 230V unit with Deep Basket



Plier Rack Stand
DDUS60P



Remote Timer
DUC315115



Drain Filter (3)
DDUS60F



Deep Basket
DDUS60DB

BioSonic® UC125 Ultrasonic Cleaning System



**UC125 Countertop Unit
Outside Dimensions**
16.5" W x 10.5" D x 13.75" H
(41.9 cm x 26.7 cm x 34.9 cm)
Height with lid open - 23-1/2"

Tank Dimensions
11.25" W x 5.625"D x 5"H
(28.6cm x 14.2cm x 12.7cm)

Tank Capacity
1.25 Gallons

Operating Frequency
55 kHz

Shipping Weight
20 lbs

Power Output
100 Watts

UC125 unit include:
Hinged Cover
15-minute timer
Front drain control
Instrument Basket
Two year Warranty



DMP-USP

MPUS-Plus Ultrasonic Cleaning Solution

MPUS-Plus Ultrasonic Cleaning Solution is specifically formulated for DryHeat processing. It does not require rinsing prior to loading in the sterilizer and evaporates cleanly with virtually no residue or odor. An effective rust inhibitor protects your instruments. Concentrated solution dilutes 1/2 oz. per gallon of tap water. Carton of four 8-oz bottles makes 64 gallons of use-strength solution.
(Carton of 4)




DUC39

Ultrasonic Cement Remover Concentrate

Cleans temporary and permanent cements such as glass ionomer and zinc phosphate quickly and easily.

16 oz bottle makes ??? pints



MADE IN USA

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

25

DDS7000 Orthodontic Dry Heat Sterilizer

FDA 510(k) cleared as a medical sterilization device.
Completely inert process does not damage hinged instruments,
sharp edges, or non-stainless components.
Sterilizing temperature of 375°F (191°C).
Cool down cycle allows immediate handling after sterilization.
Complete cool-to-cool cycle.
Large load capacity.
Plier loading racks included in FREE accessory kit.

Capacity
36 pliers

Outside Dimensions
18-3/4" W x 19-3/4" D x 23-1/8" H
(47.6 cm x 50.1 cm x 58.7 cm)

Weight
90 lbs

Chamber Dimensions
12-1/2" W x 9" D x 6" H
(31.7 cm x 22.8 cm x 15.2 cm)

Power Requirement
15A/115V or 10A/230V



115 V DDS7000115
230 V DDS7000230

4

Sterilize

for orthodontic applications, dry heat is the preferred method of choice.


"Heat-tolerant dental instruments usually are sterilized by 1) steam under pressure (autoclaving), 2) dry heat, or 3) unsaturated chemical vapor. All sterilization should be performed by using medical sterilization equipment cleared by FDA. The sterilization times, temperatures, and other operating parameters recommended by the manufacturer of the equipment used, as well as instructions for correct use of containers, wraps, and chemical or biological indicators, should always be followed." CDC Guidelines

Accessory Kit

D7107



Includes 2 D7000C Racks, 2 D7000 Silver Racks,
1 Spore Strip Holder, 8 oz bottle of MPUS-Plus,
Sample Barrier Bag, 1 Spore Test, 1 DSY-20 Silicon Lubricant




Vent Adaptors

DAD-BL

Straight Back - mounts to sterilizer

DAD-UL

Universal - mounts to exhaust stack




Vertical Plier Rack 9 Instruments

D7000-Silver




Vertical Plier Rack Mini

D7000M-S




Horizontal Tray

D5000S-T




Standard Combo Rack

D7000C



Long Handle Combo Rack

D7000CL



Horizontal Half Tray

D5000S-HT

DSL-16
16oz Spray



DSY-20
5.5g syringe



Pure Silicon Lubricant

STERILIZER ACCESSORIES

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

Barrier Bags



Large
(fits large racks-250/pkg)

Size
DBB-101 10 1/2" X 9"

Medium
(fits mini racks-250/pkg)

Size
DBB-102 7" X 7 1/2"

Small
(fits scalars, mouth mirrors-150/pkg)

Size
DBB-103 3 1/2" X 8 1/2"

5

Store

instruments can be organized by type or procedure.
Use barrier bags to prevent airborne contamination.

"Critical and semi-critical instruments or containers of critical and semi-critical instruments sterilized by a heat or vapor method shall be packaged or wrapped before sterilization if they are not to be used immediately after being sterilized. These packages or containers shall remain sealed unless the instruments within them are placed onto a setup tray and covered with a moisture impervious barrier on the day the instruments will be used and shall be stored in a manner so as to prevent contamination." California Board of Dental Examiners Guidelines

Color Code Rings

DCCR
DCCR-K



Sterilizer Monitoring

"Biological indicators (BIs) (i.e., spore tests) are the most accepted method for monitoring the sterilization process because they assess it directly by killing known highly resistant microorganisms, rather than merely testing the physical and chemical conditions necessary for sterilization... Correct functioning of sterilization cycles should be verified for each sterilizer by the periodic use (at least weekly) of BIs." CDC Guidelines

Passport
Postage Paid Mailer
2 Strips Plus Control



DPS012

DPS052

EMS
Self Mailer
1 Strip Plus Control



DEMS012

DEMS052

SporeView
In Office Culturing



DCS020




Incubator, Log Book,
Stand, Strips

37° Starter Kit
DSK037

56° Starter Kit
DSK056



Spore Strip Holder
DSSH-1



Dry Heat Process Indicators
DIS100

STERILIZER ACCESSORIES

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com



Ultracare Disinfectant Cleaner Concentrate

EPA-registered hospital grade disinfectant.
Fungicidal. Germicidal. Pseudomonacidal. Virucidal.
One step cleaning and disinfecting.
Effective against 141 organisms in 10-minute contact time including HBV, HCV, and SARS.
Formulated with a premium surfactant.
Effective in the presence of 98% organic soil.
Complies with CDC, ADA, and OSHA guidelines.
Will not damage plastics or stainless steels.
Active ingredients Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride. Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride.

16 oz makes 32 quarts!

Spray Kit

Includes 16 oz concentrate and pint spray bottle (makes 64 pints)

DUL55SP

Starter Kit

Includes 2 oz concentrate and pint spray bottle (makes 8 pints)

DUL55SK

16 oz concentrate

bottle makes 32 quarts

DUL55

64 oz concentrate

bottle makes 128 quarts

DUL555

Ultra-Tray Disinfectant Soak Tray

1 gallon Capacity
13-1/2" x 8-3/4" x 4-3/4"

Includes self-draining basket and lid. Safe for use with all disinfectants including phenols.

DST1



MetriCide® Glutaraldehyde Disinfectant /Sterilant



DMX2800

Long-life (28-day) high level disinfectant / sterilant.
Activated Alkaline formula (Ph approximately 8)
Virucidal and bacteriocidal in 10 minutes.
Tuberculocidal in 90 minutes.
Sterilant in 10 hours.
Pleasant pine fragrance
Should not be used with orthodontic pliers.

1 Gal.



Metri-Test Strips
Measures Minimum Effective
Concentration (MEC) of
Glutaraldehydes

DMX204

10-minute Instrument Disinfectant Concentrate



D70045

Phenol-based Intermediate level disinfectant.
Immersion solution only
Effective rust inhibitor
1:128 dilution
Tuberculocidal in 10 minutes.
Do not use on rubber and plastics.

1 Pint (makes 64 quarts)

to order
1.800.523.5944
fax. 330.916.7333

info@dentronix.com
dentronix.com

Placing an Order

Phone

Toll Free 800-523-5944 USA
(330) 916-7300 Worldwide
8:00AM to 5:00PM EST to speak
with a knowledgeable Technical
Sales Representative.

FAX

24 Hours a Day
FAX Line (330) 916-7333

E-mail

E-mail us directly 24 hours a day at:
info@dentronix.com (export)
sales@dentronix.com (USA)

Internet

Visit us and shop online at:
www.dentronix.com

100% Customer satisfaction

Dentronix offers a full line of fine
instruments and equipment for the
orthodontic profession. Known for
our quality and service, Dentronix
stands behind every product with
our 100% customer satisfaction
guarantee.

Student/University Large Group Volume Purchases

Dentronix is pleased to offer special
discounts for:

- High Volume
- Students, recent graduates and
faculty members
- Study clubs and group purchases.
Please call your Sales Representa-
tive today to learn more about
these discounts.

Shipping World-wide

Dentronix products are shipped
daily throughout the world from our
manufacturing facility in Ohio.
To facilitate shipments to many of
our International customers, we
have established a network of dis-
tributors in countries outside of the
US. To find out if you are serviced by
a distributor, contact our Interna-
tional Sales Department today, or
email: info@dentronix.com.

Payment Terms

We accept:
AMEX VISA MASTERCARD

All shipments are F.O.B. Cuyahoga
Falls, OH, unless otherwise stated in
writing. Payment of invoice can be
made via check, MC, Visa, American
Express or wire transfer. (Please do
not send cash through the mail).
A 1-1/2% finance charge per month
will be applied to all accounts that
are past due the NET 30 day terms.
Note: Minimum order is \$25 or a
\$5 handling fee is applied.

Returns & Credits

Returns must be received within 30
days of invoice in restockable condi-
tion for full credit. Custom or modi-
fied manufactured products cannot
be returned for credit or exchange.
Please make sure to insure all pack-
ages being returned to Dentronix
for full value. No cash refunds after
30 days of invoice date. After 30
days a 20% restocking fee is ap-
plied to all product returns that can
be resold as new.
Prior approval must be obtained
from Dentronix before product is
returned for Credit.

Product returned without a return
authorization number - (RA #) will
be returned to the customer at their
expense. The RA # must be on the
outside of the shipping package.
Shipment of the returned product
must be prepaid.

Repairs

Dentronix will sharpen or reinsert
our instruments at a nominal
charge. Repair services are provided
on Dentronix manufactured instru-
ments only. All other instruments
will be returned to the customer.
Inserted plier repairs are warranted
for one year. Non-inserted instru-
ment repairs are warranted for 90
days. Allow 3 to 5 weeks (plus ship-
ping time) for repairs. Pliers must
be sterilized before they are sent in
for repair.
Please send your insured package
to:

Repair Department
Dentronix Inc.
235 Ascot Parkway
Cuyahoga Falls, OH 44223

Note: Repair Forms are now
available on our web site www.
dentronix.com

Plier Repair

- Sharpen
- Insert Replacement and
Reconditioning
- General Cleaning, Polishing
and Sandblasting
- Warranty repair

Adams Plier	15	Compound Action Cutter	15	How Plier, Angled	10	Plier Rack, Ultrasonic	24	Sterilizer Accessory Kit	26
Anterior Band Remover	12	Consultation Models	20	Impression Trays	20	Plier Racks	21	Sterilizer Monitors	29
Arch Contouring Plier, Non-grooved	9	Cotton Plier	19	Inserting Plier	10	Pluggers	16	Sterilizer Plier Rack	27
Arch Contouring Plier, Grooved	9	Crimpable Hook Plier	11	Jarabak Plier	5	Positioner, Bracket	13	Sterilizer Tray	27
Arch Detailing Plier	9	Crown and Collar Scissors	18	Laboratory Pliers	14, 15	Posterior Band Remover	12	Sterilization Protocol	22
Arch Plier, Tweed	7	Cutter, Compound Action	15	Ligature Cutter	4	Process Indicators	29	Stop Plier, V-Bend	8
B.I.'s	29	Cutter, Distal End	3	Light Wire and Arch Contouring Plier	7	Protocol, Sterilization	22	Test Strips, Glutaraldehyde	31
Bags, Barrier	28	Cutter, Hard Wire	4	Light Wire Plier	6	Protractor, Cephalometric	20	Three Prong Plier	14, 15
Band Contouring Plier	13	Cutter, Lab	15	Light Wire Plier w/ Cutter	6	Pusher, Band	13	Three Prong Plier, NiTi	8
Band Driver Tip	16	Cutter, Ligature	4	Light Wire, Magness	5	Pushers	16	Timer, Ultrasonic	24
Band Driver, Roth	16	DDS7000	26	Light Wire, Roth	7	Quad Helix Plier	11	Tips, Band Remover	12
Band Pusher	13	Debonding Plier	12	Lingual Arch Bending Plier	8	Racks, Sterilizer	27	Tips, Bracket Remover	12
Band Remover	12	Deep Basket, Ultrasonic	24	Lingual Arch Placing Plier	11	Racks, Plier	21	Torquing Set	8
Band Slitter	12	Detailing Plier	9	Lingual Arch Placing Plier	11	RES Expansion Screw	15	Tray, Disinfectant Soak	30
Barrier Bags	28	Dickenson Plier	13	Lingual Utility Plier	11	Ribbon Arch Plier	7	Tray, Impression	20
Bernard Plier	15	Disinfectant Products	30, 31	Loop Impingement Plier	9	Rings, Color Code	28	Tray, Sterilizer	27
Biological Monitors	29	Disinfectant, High-level	31	Loop Plier, Nance	7	Roth Band Driver	16	Turbyfill Plier	9
Bird Beak Plier	5, 14	Distal End Cutter	3	Lubricant, Plier	27	Roth Light Wire	7	Tweed Arch Plier	7
Bird Beak Plier, Serrated	5, 14	Drain Filter, DDUS60	24	Magness Light Wire	5	Safety Hold Cutter	3	Tweed Omega Loop Plier	7
Bird Beak w/ Cutter	5	Dry Heat Sterilizer	26	Mathieu Plier	17	Safety Key, RES	15	Tweezer, Bracket Placing	13
Bird Beak, Large	14	Elastic Placers	17, 19	Metricide	31	Samson Plier	7	Ultracare Disinfectant	30
Bite Stick	13	Expansion Screw, Ratchet (RES)	15	Mini-rack, Sterilizer	27	Scalers	16	Ultrasonic Cleaning Machine	23-25
Boone Positioning Gauge	13	Explorers	16	Mirror Handles	18	Scissors	18	Ultrasonic Cleaning Machine, Recessed	24
Bowls, Mixing	21	Facebow Pliers	15	Mirrors	18	Scissors, Crown and Collar	18	Ultrasonic Cleaning Solution	25
Bracket Placer	13	Festooning Scissors	18	Mixing Bowls	21	Separating Plier	19	Ultrasonic Deep Basket	24
Bracket Removing Plier	12	Glutaraldehyde	31	Models, Consultation	20	Serrated Bird Beak	5, 14	Ultrasonic Plier Rack/ Stand	24
Cam-Lok Elastic Positioner	19	Glutareldehyde Test Strips	31	Mosquito Hemostat	17	Shear Hold	3	Utility Plier	10, 11
Cement Remover Concentrate	25	Hard Wire Cutter	4	MPUS Cleaning Solution	25	Silicon Lubricant	27	V Bend Plier	8
Cephalometric Protractor	20	Height Gauge	13	Nance Loop Plier	7	Slitter, Band	12	Vent Adaptor	27
Cheek Retractors	19	Hemostats	17	NiTi Three Prong Plier	8	Soak Tray	30	Weingart Plier	10
Cinch Back Plier	11	Holder, Spore Strip	29	Omega Loop Plier, Tweed	7	Spore Strip Holder	29	Young Plier	14
Cleaner, Ultrasonic	23-25	Hook Crimping Plier	11	Optical Plier	8	Spore Strips	29		
Color Code Rings	28	How Plier	10	Parallel Action Plier	15	Step Plier	9		
				Pin Inserting Plier	10				



Instrumentation. Sterilization. Innovation.

235 Ascot Parkway

Cuyahoga Falls, OH 44223

1.800.523.5944