

Lewanit - Dental Wires



Lewanit Dental Wires are made of high quality Chrome-Nickel-Alloy. The bright shining wires are available in different qualities and quantities.

Steadyness:

Quality "soft"	(L/V 7..)	600 - 800 N/mm ²
Quality "hard"	(L/V 8..)	1400 - 1600 N/mm ²
Quality "spring hard"	(L/V 9..)	1800 - 2000 N/mm ²

Lewanit - Dental Wires

Laboratory coil

Durchmesser Diameter		Inhalt content	Artikel Nr. Article Number		
mm	inch	Meter	weich soft	hart hard	federhart spring hard
0,40	.016	50	L 704	L 804	L 904
0,50	.020	50	L 705	L 805	L 905
0,60	.024	40	L 706	L 806	L 906
0,70	.028	30	L 707	L 807	L 907
0,80	.032	20	L 708	L 808	L 908
0,90	.036	10	L 709	L 809	L 909
1,00	.040	10	L 710	L 810	L 910
1,10	.044	10	L 711	L 811	L 911
1,20	.048	10	L 712	L 812	L 912

Stock coil

Durchmesser Diameter		Inhalt content	Artikel Nr. Article Number		
mm	inch	kg	weich soft	hart hard	federhart spring hard
0,40	.016	0,5	V 704	V 804	V 904
0,50	.020	0,5	V 705	V 805	V 905
0,60	.024	0,5	V 706	V 806	V 906
0,70	.028	0,5	V 707	V 807	V 907
0,80	.032	0,5	V 708	V 808	V 908
0,90	.036	0,5	V 709	V 809	V 909
1,00	.040	0,5	V 710	V 810	V 910
1,10	.044	0,5	V 711	V 811	V 911
1,20	.048	0,5	V 712	V 812	V 912