



d'art™

Mini Ceramic Bracket




## Mini Ceramic bracket system That Combines.....


**MiniSize**  
Smaller than most of Ceramic brackets in the market

**Dove Tail Grooves**  
for pure mechanical bonding

**Low Profile Design**  
for comfort and minimum occlusal interferences



d'art™  
Mini Ceramic Bracket



E-mail : modortho@gmail.com

Website : www.modortho.com

d'art™

**Roth Prescription**

Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+12°	+5°	.124"	R	BTD-R8-001	BTD-R2-001
				L	BTD-R8-002	BTD-R2-002
2 2 Lateral	+8°	+9°	.124"	R	BTD-R8-003	BTD-R2-003
				L	BTD-R8-004	BTD-R2-004
3 3 Cuspid with hook	0°	+11°	.124"	R	BTD-R8-005	BTD-R2-005
				L	BTD-R8-006	BTD-R2-006
3 3 Cuspid	0°	+11°	.124"	R	BTD-R8-007	BTD-R2-007
				L	BTD-R8-008	BTD-R2-008
4 4 Bicuspid	-7°	0°	.124"	U	BTD-R8-009	BTD-R2-009
				L	BTD-R8-010	BTD-R2-010
5 5 Bicuspid with hook	-7°	0°	.124"	L	BTD-R8-011	BTD-R2-011

**Patient Kits**

Description	018" Slot	.022" Slot
Upper/Lower 5x5 No Hooks	BKD-88-100	BKD-82-100
Upper/Lower 5x5 Hooks on Cuspids	BKD-88-102	BKD-82-102
Upper/Lower 5x5 Hooks 3, 4 & 5	BKD-88-104	BKD-82-104
Upper/Lower 3x3 Hooks, 3	BKD-88-106	BKD-82-106
Upper/Lower 3x3 No Hooks	BKD-88-108	BKD-82-108

**MBT Prescription**

Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	+17°	+4°	.124"	R	BTD-M8-101	BTD-M2-101
				L	BTD-M8-102	BTD-M2-102
2 2 Lateral	+10°	+8°	.124"	R	BTD-M8-103	BTD-M2-103
				L	BTD-M8-104	BTD-M2-104
3 3 Cuspid with hook	0°	+9°	.124"	R	BTD-M8-105	BTD-M2-105
				L	BTD-M8-106	BTD-M2-106
3 3 Cuspid	0°	+9°	.124"	R	BTD-M8-107	BTD-M2-107
				L	BTD-M8-108	BTD-M2-108
3 3 Cuspid with hook	-7°	+9°	.124"	R	BTD-M8-109	BTD-M2-109
				L	BTD-M8-110	BTD-M2-110
4 4 Cuspid	-7°	+9°	.124"	R	BTD-M8-111	BTD-M2-111
				L	BTD-M8-112	BTD-M2-112
4 4 Bicuspid with hook	-7°	0°	.124"	R	BTD-M8-113	BTD-M2-113
				L	BTD-M8-114	BTD-M2-114
5 5 Bicuspid	-7°	0°	.124"	U	BTD-M8-115	BTD-M2-115

**Patient Kits**

Description	018" Slot	.022" Slot
Upper/Lower 5x5 No Hooks	BKD-88-100	BKD-82-100
Upper/Lower 5x5 Hooks on Cuspids	BKD-88-102	BKD-82-102
Upper/Lower 5x5 Hooks 3, 4 & 5	BKD-88-104	BKD-82-104
Upper/Lower 3x3 Hooks, 3	BKD-88-106	BKD-82-106
Upper/Lower 3x3 No Hooks	BKD-88-108	BKD-82-108

**Standard Edgewise Prescription**

Maxillary Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
1 1 Central	0°	0°	.124"	R	BTD-E8-401	BTD-E2-401
				L	BTD-E8-402	BTD-E2-402
2 2 Lateral	0°	0°	.124"	R	BTD-E8-403	BTD-E2-403
				L	BTD-E8-404	BTD-E2-404
3 3 Cuspid	0°	0°	.124"	R	BTD-E8-405	BTD-E2-405
				L	BTD-E8-406	BTD-E2-406
3 3 Cuspid with hook	0°	0°	.124"	R	BTD-E8-407	BTD-E2-407
				L	BTD-E8-408	BTD-E2-408
4 4 Bicuspid	0°	0°	.124"	U	BTD-E8-409	BTD-E2-409
				L	BTD-E8-410	BTD-E2-410
5 5 Bicuspid with hook	0°	0°	.124"	R	BTD-E8-411	BTD-E2-411
				L	BTD-E8-412	BTD-E2-412

**Patient Kits**

Description	018" Slot	.022" Slot
Upper/Lower 5x5 No Hooks	BKD-88-100	BKD-82-100
Upper/Lower 5x5 Hooks on Cuspids	BKD-88-102	BKD-82-102
Upper/Lower 5x5 Hooks 3, 4 & 5	BKD-88-104	BKD-82-104
Upper/Lower 3x3 Hooks, 3	BKD-88-106	BKD-82-106
Upper/Lower 3x3 No Hooks	BKD-88-108	BKD-82-108

**Roth Prescription**

Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
2 2 Lower Anterior	0°	0°	.119"	U	BTD-M8-119	BTD-M2-119
				L	BTD-M8-120	BTD-M2-120
3 3 Cuspid with hook	-11°	+5°	.124"	R	BTD-M8-121	BTD-M2-121
				L	BTD-M8-122	BTD-M2-122
3 3 Cuspid	-11°	+5°	.124"	R	BTD-M8-123	BTD-M2-123
				L	BTD-M8-124	BTD-M2-124
4 4 1" Bicuspid	-17°	0°	.124"	R	BTD-M8-125	BTD-M2-125
				L	BTD-M8-126	BTD-M2-126
4 4 1" Bicuspid with hook	-17°	0°	.124"	R	BTD-M8-127	BTD-M2-127
				L	BTD-M8-128	BTD-M2-128
5 5 2" Bicuspid	-22°	0°	.124"	R	BTD-M8-129	BTD-M2-129
				L	BTD-M8-130	BTD-M2-130
5 5 2" Bicuspid with hook	-22°	0°	.124"	R	BTD-M8-131	BTD-M2-131
				L	BTD-M8-132	BTD-M2-132

\*The Modern Orthodontics version of the Roth is not implied to be an exact version of any other system, nor do we claim any endorsement of Dr. Roth.

**MBT Prescription**

Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
2 2 Lower Anterior	-6°	0°	.119"	U	BTD-M8-119	BTD-M2-119
				L	BTD-M8-120	BTD-M2-120
3 3 Cuspid with hook	-6°	+3°	.124"	R	BTD-M8-121	BTD-M2-121
				L	BTD-M8-122	BTD-M2-122
3 3 Cuspid	-6°	+3°	.124"	R	BTD-M8-123	BTD-M2-123
				L	BTD-M8-124	BTD-M2-124
4 4 1" Bicuspid	-12°	+2°	.124"	R	BTD-M8-125	BTD-M2-125
				L	BTD-M8-126	BTD-M2-126
4 4 1" Bicuspid with hook	-12°	+2°	.124"	R	BTD-M8-127	BTD-M2-127
				L	BTD-M8-128	BTD-M2-128
5 5 2" Bicuspid	-17°	+2°	.124"	R	BTD-M8-129	BTD-M2-129
				L	BTD-M8-130	BTD-M2-130
5 5 2" Bicuspid with hook	-17°	+2°	.124"	R	BTD-M8-131	BTD-M2-131
				L	BTD-M8-132	BTD-M2-132

\*The Modern Orthodontics version of the McLaughlin, Bennett, and Trevis system is not implied to be an exact version of any other system, nor do we claim any endorsement of Doctors McLaughlin, Bennett, and Trevis.

**Standard Edgewise Prescription**

Mandibular Arch	Torq.	Ang.	Width	R/L	.018" Slot	.022" Slot
2 2 Lower Anterior	0°	0°	.119"	U	BTD-E8-409	BTD-E2-409
				L	BTD-E8-410	BTD-E2-410
3 3 Cuspid	0°	0°	.124"	R	BTD-E8-411	BTD-E2-411
				L	BTD-E8-412	BTD-E2-412
3 3 Cuspid with hook	0°	0°	.124"	R	BTD-E8-413	BTD-E2-413
				L	BTD-E8-414	BTD-E2-414
4 4 Bicuspid	0°	0°	.124"	U	BTD-E8-415	BTD-E2-415
				L	BTD-E8-416	BTD-E2-416
5 5 Bicuspid with hook	0°	0°	.124"	R	BTD-E8-417	BTD-E2-417
				L	BTD-E8-418	BTD-E2-418



Bonds & De-Bonds Just Like A Metal Bracket



De-bonding tool to be used



Wrong method of De-bonding



Grasp the brackets on the mesio-distal side close to the tooth



Twist the bracket 90° clockwise to safely remove the bracket